
MATERIAL DATA SHEET

Poly-SUMO-2 chains (2-8)**Cat. # ULC-210**

Poly-SUMO-2 chains can be used to investigate mechanisms of chain recognition, binding and hydrolysis by SUMO-specific isopeptidases (SENPs), SUMO-specific E3 ligases or other proteins that contain SUMO-2 binding domains. This product is formed enzymatically with wild-type human recombinant SUMO-2 linked via lysine 11 which is the point of attachment for the C-terminal glycine of the preceding SUMO-2. Mono-SUMO-2 has been removed from the chain mixture.

Product Information

Quantity:	25 µg
Stock:	X mg/ml in 50 mM Hepes pH 8.0, 100 mM NaCl, 1 mM DTT. Concentration will vary with specific Lot #.
Purity:	> 95 % by SDS-PAGE

Use & Storage

Use:	Typical concentrations will depend on specific assay conditions and method of detection.
Storage:	Store at -80°C. Avoid multiple freeze/thaw cycles.

Literature

References:	Bylebyl G.R., <i>et al.</i> (2003) <u>J. Biol. Chem.</u> 278 :44113-44120 Dohmen R. J. (2004) <u>Biochem. Biophys. Acta.</u> 1695 :113-131 Johnson E. (2004) <u>Ann. Rev. Biochem.</u> 73 :355-382 Tatham M.H., <i>et al.</i> (2001) <u>J. Biol. Chem.</u> 276 : 35368-35374
--------------------	---

For Laboratory Research Use Only, Not For Use in Humans