

MATERIAL DATA SHEET

Recombinant Human Ubiquitin Mutant with K29 only

Cat. # UM-K290

Ubiquitin is a 76 amino acid (aa) protein that is ubiquitously expressed in all eukaryotic organisms. Ubiquitin is highly conserved with 96% aa sequence identity shared between human and yeast Ubiquitin, and 100% aa sequence identity shared between human and mouse Ubiquitin (1). In mammals, four Ubiquitin genes encode for two Ubiquitin-ribosomal fusion proteins and two poly-Ubiquitin proteins. Cleavage of the Ubiquitin precursors by deubiquitinating enzymes gives rise to identical Ubiquitin monomers each with a predicted molecular weight of 8.6 kDa. Conjugation of Ubiquitin to target proteins involves the formation of an isopeptide bond between the C-terminal glycine residue of Ubiquitin and a lysine residue in the target protein. This process of conjugation, referred to as ubiquitination or ubiquitylation, is a multi-step process that requires three enzymes: a Ubiquitin-activating (E1) enzyme, a Ubiquitin-conjugating (E2) enzyme, and a Ubiquitin ligase (E3). Ubiquitination is classically recognized as a mechanism to target proteins for degradation and as a result, Ubiquitin was originally named ATP-dependent Proteolysis Factor 1 (APF-1) (2,3). In addition to protein degradation, ubiquitination has been shown to mediate a variety of biological processes such as signal transduction, endocytosis, and post-endocytic sorting (4-7).

This Ubiquitin mutant contains only a single lysine, K29, with all other lysines mutated to arginine. This mutation renders Ubiquitin able to form poly-Ubiquitin isopeptide chains with other Ubiquitin molecules only via the K29 lysine.

Product Information

Quantity:	1 mg
MW:	8.7 kDa
Source:	<i>E. coli</i> -derived Accession # P0CG47
Stock:	Lyophilized from a solution in deionized water.
Solubility:	Reconstitute at 10 mg/mL in an aqueous solution.
Purity:	>95%, by SDS-PAGE under reducing conditions and visualized by Colloidal Coomassie® Blue stain.

Use & Storage

- Use:** The lysine residue utilized for Ubiquitin chain formation is functionally important. Single lysine Ubiquitin mutants are ideal for investigating biological processes involving specific Ubiquitin chain linkage. Recombinant Human Ubiquitin K29-only allows for the formation of only K29-linked Ubiquitin chains. Reaction conditions will need to be optimized for each specific application. We recommend an initial Recombinant Human Ubiquitin K29-only concentration of 0.2-1 mM.
- Storage:** **Use a manual defrost freezer and avoid repeated freeze-thaw cycles.**
- 12 months from date of receipt, -20 to -70 °C as supplied.
 - 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Literature

References:

1. Sharp, P.M. & W.-H. Li. (1987) Trends Ecol. Evol. **2**:328.
2. Ciechanover, A. *et al.* (1980) Proc. Natl. Acad. Sci. USA **77**:1365.
3. Herskho, A. *et al.* (1980) Proc. Natl. Acad. Sci. USA **77**:1783.
4. Greene, W. *et al.* (2012) PLoS Pathog. **8**:e1002703.
5. Tong, X. *et al.* (2012) J. Biol. Chem. **287**:25280.
6. Wei, W. *et al.* (2004) Nature **428**:194.
7. Wertz, I.E. *et al.* (2004) Nature **430**:694.
8. Arnason T., *et al.* (1994) Mol. Cell. Biol. **14**:7876-7883.
9. Spence J., *et al.* (1995) Mol. Cell. Biol. **15**:1265-1273.
10. Russel N.S., *et al.* (2004) Biochemistry **43**:4844-4854.

For research use only. Not for use in humans.