



***Recombinant Mouse BLK  
Western Blotting Standard***

**Catalog Number:** WBC007

**Size:** 100  $\mu$ L

**Standard Preparation:** Crude extracts of *E. coli*-expressing amino acids 1 - 124 of recombinant mouse BLK were solubilized in SDS sample buffer (3% SDS, 0.125 M Tris pH 6.8, 5% glycerol, and bromophenyl blue) containing 10 mM dithiothreitol.

**Western Blot:** Heat sample at 90° - 100° C for 30 - 60 seconds and use 10  $\mu$ L per lane.

**Storage:** Store at -20° C.