



***Recombinant Human Hrk
Western Blotting Standard***

Catalog Number: WBC013

Size: 100 μ L

Standard Preparation: *E. coli* expressing full-length (amino acids 1 - 91) recombinant human harakiri (Hrk), an activator of apoptosis, fused to green fluorescent protein (GFP) is solubilized in SDS sample buffer (3% SDS, 0.125 M Tris pH 6.8, 5% glycerol, and bromophenyl blue) containing 10 mM dithiothreitol. Hrk-GFP electrophoreses as a 35 kDa protein.

Western Blot: Heat sample at 90° - 100° C for 30 - 60 seconds and use 10 μ L per lane.

Storage: Store at -20° C.