



***Recombinant Human/Mouse MMP-24
Western Blotting Standard***

Catalog Number: WBC014

Size: 25 Lanes

Standard Preparation: Medium conditioned by MDCK cells expressing the extracellular domain of recombinant mouse MMP-24 (MT5-MMP) (Wang, *et al.*, 1999, FEBS Lett. **462**:261 - 266) was aliquoted and lyophilized. MMP-24 standard electrophoreses as a 58 kDa protein.

Western Blot: Resuspend the lyophilized powder in 250 μ L SDS-PAGE sample buffer. Load 10 μ L per lane.

Storage: Store at -20° C.