



***Recombinant Human TIMP-2  
Western Blotting Standard***

**Catalog Number:** WBC023

**Size:** 25 Lanes

**Standard Preparation:** Conditioned media from CHO cells expressing recombinant human tissue inhibitor of metalloproteinases-2 (rhTIMP-2) was aliquoted and lyophilized. In western blots, using R&D Systems' anti-human TIMP-2 (Catalog # AF971 and MAB971), rhTIMP-2 standard can be detected as a 25 kDa protein under both reducing and non-reducing conditions.

**Western Blot:** Resuspend the lyophilized powder in 250  $\mu$ L SDS-PAGE sample buffer. Load 10  $\mu$ L per lane.

**Storage:** Store at -20° C.