



***Recombinant Human MMP-1
Western Blotting Standard***

Catalog Number: WBC024

Size: 25 Lanes

Standard Preparation: Conditioned media from NS0 mouse myeloma cells expressing recombinant human matrix metalloproteinase-1 (rhMMP-1) was aliquoted and lyophilized. In western blots, using R&D Systems' anti-human MMP-1 antibodies (Catalog # MAB900, MAB901 and AF901), rhMMP-1 standard can be detected as a 55 kDa protein under both reducing and non-reducing conditions.

Western Blot: Resuspend the lyophilized powder in 250 μ L SDS-PAGE sample buffer. Load 10 μ L per lane.

Storage: Store at -20° C.