



***Recombinant Human TACE Ectodomain
Western Blotting Standard***

Catalog Number: WBC029

Size: 25 Lanes

Standard Preparation: Conditioned media from Sf21 insect cells expressing the ectodomain of recombinant human TNF- α converting enzyme (rhTACE (ecto)) was aliquoted and lyophilized. In western blots using R&D Systems' anti-human TACE (ecto) antibody (Catalog # MAB930), rhTACE (ecto) standard can be detected as a 75 kDa protein under non-reducing conditions.

Western Blot: Resuspend the lyophilized powder in 250 μ L SDS-PAGE sample buffer. Load 10 μ L per lane.

Storage: Store at -20° C.