



Recombinant Mouse MMP-9 Western Blotting Standard

Catalog Number: WBC031

Size: 25 Lanes

Standard Preparation: Conditioned media from NS0 mouse myeloma cells expressing recombinant mouse matrix metalloproteinase-9 (rmMMP-9) was aliquoted and lyophilized. In western blots, using R&D Systems' anti-mouse MMP-9 (Catalog # BAF909), rmMMP-9 standard can be detected as a 110 kDa protein under both reducing and non-reducing conditions.

Western Blot: Resuspend the lyophilized powder in 250 μ L SDS-PAGE sample buffer. Load 10 μ L per lane.

Storage: Store at ≤ -20 °C.