

Product Datasheet

Lymphotoxin beta R/TNFRSF3 Antibody [DyLight 680] AF1008Y

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF1008Y

Updated 7/7/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF1008Y



AF1008Y

Lymphotoxin beta R/TNFRSF3 Antibody [DyLight 680]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 680
Purity	Antigen Affinity-purified
Buffer	50mM Sodium Borate
Product Description	
Host	Goat
Gene ID	4055
Gene Symbol	LTBR
Species	Mouse
Specificity/Sensitivity	Detects mouse Lymphotoxin β R/TNFRSF3 in direct ELISA and Western blots. In direct ELISAs, approximately 5% cross-reactivity with recombinant human (rh) LT β R is observed, and less than 1% cross-reactivity with recombinant mouse (rm) 4-1BB, rmCD27, rmOPG, rmRANK, rmCD30, rmCD40, rhDR6, rmTNF RI, rmTNF RII, rmEDAR, rmFas and rmGITR is observed.
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse Lymphotoxin β R/TNFRSF3 Ser28-Pro218 Accession # P50284
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF1008Y

AB-108-CY	Goat IgG Isotype Control [DyLight 680]
NBP2-51906-0.05mg	Recombinant Human Lymphotoxin beta R/TNFRSF3 His Protein
210-TA-005	TNF-alpha [Unconjugated]
7538-LR-100	Lymphotoxin beta R/TNFRSF3 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF1008Y

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

