

Product Datasheet

FGFR5/FGFRL1 Antibody [DyLight 350] AF2257D

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF2257D

Updated 7/7/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF2257D



AF2257D**FGFR5/FGFRL1 Antibody [DyLight 350]**

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Polyclonal |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | DyLight 350 |
| Purity | Antigen Affinity-purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Goat |
| Gene ID | 53834 |
| Gene Symbol | FGFRL1 |
| Species | Human |
| Specificity/Sensitivity | Detects human FGF R5/FGFRL1 in direct ELISAs and Western blots. In direct ELISAs and Western blots, approximately 65% cross-reactivity with recombinant mouse (rm) FGF R5 and rmFGF R5 β is observed and less than 1% cross-reactivity with recombinant human (rh) FGF R1 α (IIIb), rhFGF R1 α (IIIc), rhFGF R1 β (IIIb), rhFGF R1 β (IIIc), rhFGF R2 α (IIIb), rhFGF R2 α (IIIc), rhFGF R2 β (IIIb), rhFGF R2 β (IIIc), rhFGF R3 α (IIIb), rhFGF R3 (IIIc), and rhFGF R4 is observed. |
| Immunogen | Mouse myeloma cell line NS0-derived recombinant human FGF R5/FGFRL1 Ala25-Pro378 Accession # Q8N441 |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Western Blot |
| Recommended Dilutions | Western Blot |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF2257D

| | |
|---------------|--|
| AB-108-CD | Goat IgG Isotype Control [DyLight 350] |
| NBP2-56879PEP | FGFR5/FGFRL1 Recombinant Protein Antigen |
| 233-FB-025 | FGF basic/FGF2/bFGF [Unconjugated] |
| 1899-FR-050 | FGFR5/FGFRL1 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF2257D

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

