

Product Datasheet

FGF-4 Antibody [DyLight 405] AF235E

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF235E

Updated 7/7/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF235E



AF235E

FGF-4 Antibody [DyLight 405]

| Product Information | |
|-----------------------------|--|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Polyclonal |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | DyLight 405 |
| Purity | Antigen Affinity-purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Goat |
| Gene ID | 2249 |
| Gene Symbol | FGF4 |
| Species | Human |
| Specificity/Sensitivity | Detects human FGF□4 in direct ELISAs and Western blots. In direct ELISAs, approximately 20% cross-reactivity with recombinant human (rh) FGF-6 is observed and less than 2% cross-reactivity with bovine FGF acidic, bovine FGF basic, rhFGF acidic, rhFGF basic, rhβ□ECGF, rhFGF-5, rhFGF□7, rhFGF□8b, recombinant mouse (rm) FGF□8b, rmFGF□8c, and rhFGF□10 is observed. |
| Immunogen | <i>E. coli</i> -derived recombinant human FGF-4 Ala31-Leu206 Accession # P08620 |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Western Blot, Immunohistochemistry, Neutralization |
| Recommended Dilutions | Western Blot, Immunohistochemistry, Neutralization |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF235E

| | |
|----------------|--|
| AB-108-CE | Goat IgG Isotype Control [DyLight 405] |
| NBP2-34864-5ug | Recombinant Human FGF-4 Protein |
| 236-EG-200 | EGF [Unconjugated] |
| 235-F4-025 | FGF-4 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF235E

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

