

Product Datasheet

IL-18/IL-1F4 Antibody [DyLight 405] AF2924E

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF2924E

Updated 7/20/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF2924E



AF2924E

IL-18/IL-1F4 Antibody [DyLight 405]

| Product Information | |
|-----------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Polyclonal |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | DyLight 405 |
| Purity | Antigen Affinity-purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Goat |
| Gene ID | 3606 |
| Gene Symbol | IL18 |
| Species | Canine |
| Specificity/Sensitivity | Detects canine IL \square 18/IL \square 1F4 in direct ELISAs and Western blots. In direct ELISAs and Western blots, approximately 50% cross-reactivity with recombinant porcine IL-18 is observed, 20% cross-reactivity with recombinant rhesus macaque IL-18 is observed, 15% cross-reactivity with recombinant human IL-18 is observed, and 5% cross-reactivity with recombinant mouse IL-18 and recombinant rat IL-18 is observed. |
| Immunogen | <i>E. coli</i> -derived recombinant canine IL \square 18/IL \square 1F4 Tyr37-Ser193 Accession # Q9XSR0 |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Western Blot, Immunocytochemistry |
| Recommended Dilutions | Western Blot, Immunocytochemistry |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF2924E

| | |
|-----------------|--|
| AB-108-CE | Goat IgG Isotype Control [DyLight 405] |
| NBP2-59520-20ug | Recombinant Mouse IL-18/IL-1F4 Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| DL180 | IL-18/IL-1F4 [HRP] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF2924E

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

