

# Product Datasheet

## Ephrin-A5 Antibody [Janelia Fluor® 549] AF3743I

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF3743I](http://www.novusbio.com/AF3743I)

Updated 7/7/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF3743I](http://www.novusbio.com/reviews/destination/AF3743I)



**AF3743I**

Ephrin-A5 Antibody [Janelia Fluor® 549]

Product Information	
<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	Please see the vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at 4C in the dark.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	0.05% Sodium Azide
<b>Isotype</b>	IgG
<b>Conjugate</b>	Janelia Fluor 549
<b>Purity</b>	Antigen Affinity-purified
<b>Buffer</b>	50mM Sodium Borate
Product Description	
<b>Host</b>	Goat
<b>Gene ID</b>	1946
<b>Gene Symbol</b>	EFNA5
<b>Species</b>	Human
<b>Specificity/Sensitivity</b>	Detects human Ephrin-A5 in direct ELISAs and Western blots. In direct ELISAs, less than 1% cross-reactivity with recombinant mouse (rm) Ephrin-A1, rmEphrin-A2, recombinant human (rh) Ephrin-A3, rhEphrin-A4, rhEphrin-B3, rmEphrin-A4, rmEphrin-B1 and rmEphrin-B2 is observed.
<b>Immunogen</b>	Mouse myeloma cell line NS0-derived recombinant human Ephrin-A5 Gln21-Asn203 Accession # P52803
<b>Notes</b>	Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.
Product Application Details	
<b>Applications</b>	Western Blot, Immunohistochemistry
<b>Recommended Dilutions</b>	Western Blot, Immunohistochemistry
<b>Application Notes</b>	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF3743I**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
AB-108-CI	Goat IgG Isotype Control [Janelia Fluor 549]
NBP2-58251PEP	Ephrin-A5 Recombinant Protein Antigen
602-A1-200	Ephrin-A1

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF3743I](http://www.novusbio.com/reviews/submit/AF3743I)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

