

# Product Datasheet

## SNF1LK2/SIK2 Antibody [DyLight 650] AF5737W

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF5737W](http://www.novusbio.com/AF5737W)

Updated 5/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF5737W](http://www.novusbio.com/reviews/destination/AF5737W)



**AF5737W**

SNF1LK2/SIK2 Antibody [DyLight 650]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Polyclonal  |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG   |
| Conjugate           | DyLight 650   |
| Purity              | Antigen Affinity-purified   |
| Buffer              | 50mM Sodium Borate  |

| Product Description     |  |
|-------------------------|--|
| Host                    | Goat   |
| Gene ID                 | 23235  |
| Gene Symbol             | SIK2   |
| Species                 | Human  |
| Specificity/Sensitivity | Detects human SNF1LK2 in Western blots.  |
| Immunogen               | <i>E. coli</i> -derived recombinant human SNF1LK2<br>Ala663-Leu818<br>Accession # Q9H0K1 |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.        |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, Immunocytochemistry                                      |
| Recommended Dilutions       | Western Blot, Immunocytochemistry                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF5737W**

---

|                    |   |
|--------------------|---|
| AB-108-CW          | Goat IgG Isotype Control [DyLight 650]              |
| H00023235-Q01-10ug | Recombinant Human SNF1LK2/SIK2 GST (N-Term) Protein |
| DYC2510-2          | CREB [p Ser133] [Biotin]                            |
| NBL1-16279         | SNF1LK2/SIK2 Overexpression Lysate                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF5737W](http://www.novusbio.com/reviews/submit/AF5737W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

