

Product Datasheet

S100A16 Antibody [FITC] AF6397F

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF6397F

Updated 5/14/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF6397F



AF6397F**S100A16 Antibody [FITC]**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	FITC
Purity	Antigen Affinity-purified
Buffer	PBS
Product Description	
Host	Sheep
Gene ID	140576
Gene Symbol	S100A16
Species	Mouse
Specificity/Sensitivity	Detects mouse S100A16 in direct ELISAs and Western blots. In direct ELISAs, approximately 30% cross-reactivity with recombinant human S100A16 is observed, and less than 1% cross-reactivity with recombinant mouse S100A4, recombinant human (rh) S100A2, and rhS100A10 is observed.
Immunogen	E.coli-derived recombinant mouse S100A16 Ala2-Ser124 Accession # NP_080692
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF6397F

NBP1-96898	Sheep IgG Isotype Control [FITC]
NBP1-51085-0.1mg	Recombinant Human S100A16 His Protein
NBL1-15653	S100A16 Overexpression Lysate
AF887	AKT [p Ser473] Antibody [Unconjugated] - Pan Specific

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF6397F

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

