

# Product Datasheet

## **ST8 alpha-2,8-Sialyltransferase 8A/ST8SIA1/Ganglioside GD3 Antibody [DyLight 488] AF6716K**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF6716K](http://www.novusbio.com/AF6716K)

Updated 5/8/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF6716K](http://www.novusbio.com/reviews/destination/AF6716K)



**AF6716K**ST8  $\alpha$ -2,8-Sialyltransferase 8A/ST8SIA1/Ganglioside GD3 Antibody [DyLight 488]

| Product Information         |                                                                                                                                                    |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                             |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                        |
| Storage                     | Store at 4C in the dark.                                                                                                                           |
| Clonality                   | Polyclonal                                                                                                                                         |
| Preservative                | 0.05% Sodium Azide                                                                                                                                 |
| Isotype                     | IgG                                                                                                                                                |
| Conjugate                   | DyLight 488                                                                                                                                        |
| Purity                      | Antigen Affinity-purified                                                                                                                          |
| Buffer                      | 50mM Sodium Borate                                                                                                                                 |
| Product Description         |                                                                                                                                                    |
| Host                        | Sheep                                                                                                                                              |
| Gene ID                     | 6489                                                                                                                                               |
| Gene Symbol                 | ST8SIA1                                                                                                                                            |
| Species                     | Human                                                                                                                                              |
| Specificity/Sensitivity     | Detects human ST8 $\alpha$ -2,8-Sialyltransferase 8A/ST8SIA1 in direct ELISAs and Western blots.                                                   |
| Immunogen                   | Chinese hamster ovary cell line CHO-derived recombinant human ST8 $\alpha$ -2,8-Sialyltransferase 8A/ST8SIA1<br>Tyr49-Ser356<br>Accession # Q92185 |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                  |
| Product Application Details |                                                                                                                                                    |
| Applications                | Western Blot, Immunohistochemistry                                                                                                                 |
| Recommended Dilutions       | Western Blot, Immunohistochemistry                                                                                                                 |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                             |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF6716K**

---

|             |                                                                           |
|-------------|---------------------------------------------------------------------------|
| 210-TA-005  | TNF-alpha [Unconjugated]                                                  |
| 6716-GT-020 | ST8 alpha-2,8-Sialyltransferase 8A/ST8SIA1/Ganglioside GD3 [Unconjugated] |
| 202-IL-010  | IL-2 [Unconjugated]                                                       |
| DMP900      | MMP-9 [HRP]                                                               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF6716K](http://www.novusbio.com/reviews/submit/AF6716K)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

