

Product Datasheet

CLEC10A/CD301 Antibody (121A5.01) - Azide and BSA Free DDX0011P-100

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/DDX0011P-100

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/DDX0011P-100



DDX0011P-100

CLEC10A/CD301 Antibody (121A5.01) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	121A5.01
Preservative	No Preservative
Isotype	IgG1
Purity	Ion exchange chromatography
Buffer	Tris-NaCl (pH 8.0)

Product Description	
Description	Novus Biologicals Mouse CLEC10A/CD301 Antibody (121A5.01) - Azide and BSA Free (DDX0011P-100) is a monoclonal antibody validated for use in IHC, WB, Flow and IP. Anti-CLEC10A/CD301 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10462
Gene Symbol	CLEC10A
Species	Human
Reactivity Notes	Human
Immunogen	CD301 was used as immunogen.
Notes	This product is manufactured by Eurobio Scientific (formerly Dendritics) and distributed by Novus Biologicals.

Product Application Details	
Applications	Western Blot, Flow Cytometry, Flow (Cell Surface), Functional, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Western Blot 1:100-1:2000, Flow Cytometry 1:10-1:1000, Immunohistochemistry 1:10-1:500, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500, Functional, Flow (Cell Surface) 1:10-1:1000
Application Notes	Functional Activity note: internalization of the receptor following its binding by the Ab and when incubated at 37 degrees C (as per Valladeau et al. 2010)




Images

Immunohistochemistry: CLEC10A/CD301 Antibody (121A5.01) [DDX0011P-100] - Tonsil tissue section : cells stained with 121A5.01 (blue) are found throughout the tonsil except in follicles(f)



Flow Cytometry: CLEC10A/CD301 Antibody (121A5.01) [DDX0011P-100] - Monocyte-derived DC (GM+IL4) stained with 121A5.01+GAM FITC



Publications

Debeer S et al. Comparative histology and immunohistochemistry of porcine versus human skin. Eur J Dermatol 2013 -07-01 [PMID: 24047577]

Valladeau J et al. Immature human dendritic cells express asialoglycoprotein receptor isoforms for efficient receptor-mediated endocytosis. J Immunol 2001-11-01 [PMID: 11698450]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to DDX0011P-100

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/DDX0011P-100

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

