

# Product Datasheet

## OX40/TNFRSF4 Antibody (2553B) [HRP] FAB10545H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB10545H](http://www.novusbio.com/FAB10545H)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB10545H](http://www.novusbio.com/reviews/destination/FAB10545H)



**FAB10545H**

OX40/TNFRSF4 Antibody (2553B) [HRP]

| Product Information         |                                                                                                       |
|-----------------------------|-------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.           |
| Storage                     | Store at 4C in the dark.                                                                              |
| Clonality                   | Monoclonal                                                                                            |
| Clone                       | 2553B                                                                                                 |
| Preservative                | No Preservative                                                                                       |
| Isotype                     | IgG                                                                                                   |
| Conjugate                   | HRP                                                                                                   |
| Purity                      | 0                                                                                                     |
| Buffer                      | PBS                                                                                                   |
| Product Description         |                                                                                                       |
| Host                        | Rabbit                                                                                                |
| Gene ID                     | 7293                                                                                                  |
| Species                     | Human                                                                                                 |
| Specificity/Sensitivity     | Detects human OX40/TNFRSF4 in direct ELISAs.                                                          |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human OX40/TNFRSF4 Leu29-Ala216<br>Accession # P43489 |
| Product Application Details |                                                                                                       |
| Applications                | Flow Cytometry, CyTOF-ready                                                                           |
| Recommended Dilutions       | Flow Cytometry, CyTOF-ready                                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB10545H**

---

|                |                                  |
|----------------|----------------------------------|
| NBP2-24891H    | Rabbit IgG Isotype Control [HRP] |
| 210-TA-005     | TNF-alpha [Unconjugated]         |
| 1256-OX-050    | OX40/TNFRSF4 [Unconjugated]      |
| 6507-IL-010/CF | IL-4 [Unconjugated]              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB10545H](http://www.novusbio.com/reviews/submit/FAB10545H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

