

# Product Datasheet

## JNK1/2/3 [p Tyr185, p Thr183] Antibody (1006A) [HRP] FAB1205H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB1205H](http://www.novusbio.com/FAB1205H)

Updated 5/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB1205H](http://www.novusbio.com/reviews/destination/FAB1205H)



**FAB1205H**

JNK1/2/3 [p Tyr185, p Thr183] Antibody (1006A) [HRP]

<b>Product Information</b>	
<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	Please see the vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at 4C in the dark.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	1006A
<b>Preservative</b>	No Preservative
<b>Isotype</b>	IgG
<b>Conjugate</b>	HRP
<b>Purity</b>	Protein A or G purified from cell culture supernatant
<b>Buffer</b>	PBS
<b>Product Description</b>	
<b>Host</b>	Rabbit
<b>Gene ID</b>	5599
<b>Gene Symbol</b>	MAPK8IP3
<b>Species</b>	Human, Mouse, Rat
<b>Specificity/Sensitivity</b>	Detects human, mouse and rat p46 and p54 JNK when dually phosphorylated at sites homologous to T183/Y185 of JNK1 and JNK2, and T221/Y223 of JNK3 in Western blots.
<b>Immunogen</b>	Phosphopeptide containing human, rat, and mouse JNK1 T183/Y185 site
<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, Immunocytochemistry
<b>Recommended Dilutions</b>	Western Blot, Immunocytochemistry
<b>Application Notes</b>	Optimal dilution of this antibody should be experimentally determined.



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB1205H**

---

NBP2-24891H	Rabbit IgG Isotype Control [HRP]
210-TA-005	TNF-alpha [Unconjugated]
DYC1387B-2	JNK1/2/3 [Biotin]
AF467	EphB2 Antibody [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB1205H](http://www.novusbio.com/reviews/submit/FAB1205H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

