

Product Datasheet

CXCR3 Antibody (49801R) [PE] FAB160RP

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB160RP

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB160RP



FAB160RP**CXCR3 Antibody (49801R) [PE]**


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	49801R
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	PE
Purity	Protein A or G purified from cell culture supernatant
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	2833
Species	Human
Specificity/Sensitivity	Detects human CXCR3 in direct ELISAs.
Immunogen	NS0 mouse myeloma cell line transfected with human CXCR3 Met1-Leu368 Accession # P49682

Product Application Details	
Applications	Flow Cytometry, Immunohistochemistry, CyTOF-ready, Neutralization
Recommended Dilutions	Flow Cytometry, Immunohistochemistry, Neutralization, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

CXCR3 Antibody (49801R) [PE] [FAB160RP] - Vial of PE conjugated antibody. PE has two excitation maxima, 498 nm excited by the Blue laser (488 nm) and 565 nm excited by the Yellow-Green laser (561 nm). Both result in emission at 578 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB160RP

NBP1-97005PE	Mouse IgG1 Isotype Control (MG1) [PE]
NB100-56404PEP	CXCR3 Antibody Blocking Peptide
210-TA-005	TNF-alpha [Unconjugated]
H00002833-P01-2ug	Recombinant Human CXCR3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB160RP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

