

# Product Datasheet

## JNK1 Antibody (228601) [PerCP] FAB17761C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB17761C](http://www.novusbio.com/FAB17761C)

Updated 5/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB17761C](http://www.novusbio.com/reviews/destination/FAB17761C)



**FAB17761C**

JNK1 Antibody (228601) [PerCP]

| Product Information         |  |
|-----------------------------|--|
| Unit Size                   | 0.1 ml   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| Storage                     | Store at 4C in the dark.   |
| Clonality                   | Monoclonal   |
| Clone                       | 228601   |
| Preservative                | 0.05% Sodium Azide   |
| Isotype                     | IgG1   |
| Conjugate                   | PerCP  |
| Purity                      | Protein A or G purified from hybridoma culture supernatant   |
| Buffer                      | PBS  |
| Product Description         |  |
| Host                        | Mouse  |
| Gene ID                     | 5599   |
| Gene Symbol                 | MAPK8  |
| Species                     | Human, Mouse, Rat  |
| Specificity/Sensitivity     | Detects human JNK1 in direct ELISAs and human, mouse, and rat JNK1 expressed predominantly as p46 JNK (isoforms 2 and/or 3, both 384 aa), with lesser amounts of p54 JNK (isoforms 1 and/or 4, both 427 aa), in Western blots. In Western blots, no cross-reactivity with recombinant human (rh) JNK2 or rhJNK3 is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant human JNK1 isoform 2<br>Ser2-Ile384<br>Accession # P45983-1  |
| Product Application Details |  |
| Applications                | Western Blot, Immunocytochemistry  |
| Recommended Dilutions       | Western Blot, Immunocytochemistry  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.   |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB17761C](http://www.novusbio.com/reviews/submit/FAB17761C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

