

# Product Datasheet

## Histone H2AX [p Ser139] Antibody (2207D) [HRP] FAB2288H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB2288H](http://www.novusbio.com/FAB2288H)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB2288H](http://www.novusbio.com/reviews/destination/FAB2288H)



**FAB2288H**

Histone H2AX [p Ser139] Antibody (2207D) [HRP]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | 2207D                                                                                       |
| Preservative                | No Preservative                                                                             |
| Isotype                     | IgG                                                                                         |
| Conjugate                   | HRP                                                                                         |
| Purity                      | Protein A or G purified from cell culture supernatant                                       |
| Buffer                      | PBS                                                                                         |
| Product Description         |                                                                                             |
| Host                        | Rabbit                                                                                      |
| Gene ID                     | 3014                                                                                        |
| Species                     | Human                                                                                       |
| Specificity/Sensitivity     | Detects human Histone H2AX when phosphorylated at S139 in Western blots.                    |
| Immunogen                   | Phospho-peptide containing human Histone H2AX S139 site<br>Accession # P16104               |
| Product Application Details |                                                                                             |
| Applications                | Western Blot, CyTOF-ready, Intracellular Staining by Flow Cytometry                         |
| Recommended Dilutions       | Western Blot, Intracellular Staining by Flow Cytometry, CyTOF-ready                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB2288H**

---

|             |                                            |
|-------------|--------------------------------------------|
| NBP2-24891H | Rabbit IgG Isotype Control [HRP]           |
| NBP1-42464  | Recombinant Human Histone H2AX His Protein |
| DYC2288-2   | Histone H2AX [p Ser139] [Biotin]           |
| NB100-304   | 53BP1 Antibody - BSA Free                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB2288H](http://www.novusbio.com/reviews/submit/FAB2288H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

