

# Product Datasheet

## Renin Antibody (434717) [Janelia Fluor® 669] FAB4277JF669

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB4277JF669](http://www.novusbio.com/FAB4277JF669)

Updated 7/20/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB4277JF669](http://www.novusbio.com/reviews/destination/FAB4277JF669)



**FAB4277JF669**

Renin Antibody (434717) [Janelia Fluor® 669]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	434717
Preservative	0.05% Sodium Azide
Isotype	IgG2a
Conjugate	Janelia Fluor 669
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Host	Rat
Gene ID	5972
Gene Symbol	Ren1
Species	Mouse
Specificity/Sensitivity	Detects mouse Renin 1 Propeptide in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with the mature form of mouse Renin 1, recombinant human Renin, recombinant mouse (rm) Cathepsin D, or rmCathepsin E is observed.
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse Renin 1 Leu22-Arg402 Accession # P06281
Notes	Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB4277JF669**

---

NBP2-52224-0.02mg	Recombinant Mouse Renin His Protein
DAN00	Angiogenin [HRP]
DY4277	Renin [Biotin]
NBP2-12446	NLRP3/NALP3 Antibody - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB4277JF669](http://www.novusbio.com/reviews/submit/FAB4277JF669)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

