

# Product Datasheet

## TCF-3/E2A Antibody (401709) [DyLight 594] FAB44061M

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB44061M](http://www.novusbio.com/FAB44061M)

Updated 7/20/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB44061M](http://www.novusbio.com/reviews/destination/FAB44061M)



**FAB44061M**

TCF-3/E2A Antibody (401709) [DyLight 594]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	401709
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	DyLight 594
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate

Product Description	
Host	Rat
Gene ID	6929
Gene Symbol	TCF3
Species	Human
Specificity/Sensitivity	Detects human TCF-3/E2A in direct ELISAs. In this format, 100% cross-reactivity with recombinant human (rh) TCF-12 is observed and no cross-reactivity with rhTCF-1, 2, 4, 7L1, 7L2 or rhLEF-1 is observed.
Immunogen	<i>E. coli</i> -derived recombinant human TCF-3/E2A Gly100-Glu500 Accession # NP_003191
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunocytochemistry
Recommended Dilutions	Immunohistochemistry-Paraffin, Immunocytochemistry
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB44061M**

---

NBP2-38961PEP	TCF-3/E2A Recombinant Protein Antigen
NBL1-16774	TCF-3/E2A Overexpression Lysate
AF1329	beta-Catenin Antibody [Unconjugated]
NB100-56511	MyoD Antibody (5.8A) - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB44061M](http://www.novusbio.com/reviews/submit/FAB44061M)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

