

Product Datasheet

CCL21/6Ckine Antibody (59106) [Janelia Fluor® 549] FAB457I

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB457I

Updated 10/7/2024 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB457I



FAB457I

CCL21/6Ckine Antibody (59106) [Janelia Fluor® 549]

| Product Information | |
|-----------------------------|--|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 59106 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2b |
| Conjugate | Janelia Fluor 549 |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Rat |
| Gene ID | 6366 |
| Gene Symbol | Ccl21a |
| Species | Mouse |
| Specificity/Sensitivity | Detects mouse CCL21/6Ckine in direct ELISAs and Western blots. Does not cross-react with recombinant human (rh) 6Ckine, rhMIP-1 α , rhMIP-1 β , recombinant mouse (rm) MIP-1 α , rmMIP-1 β , rhMIP-1 δ , rmMIP-1 γ , rmMIP-2, rhMIP-3 α , rmMIP-3 α , rhMIP-3 β , rmMIP-3 β , rhTeck, or rmTeck. |
| Immunogen | <i>E. coli</i> -derived recombinant mouse CCL21/6Ckine Ser24-Gly133 Accession # P84444 |
| Notes | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus. |
| Product Application Details | |
| Applications | Western Blot, Immunohistochemistry, CyTOF-ready, ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair), Intracellular Staining by Flow Cytometry |
| Recommended Dilutions | Western Blot, Immunohistochemistry, ELISA Standard (Matched Pair), ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), Intracellular Staining by Flow Cytometry, CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB457I

| | |
|----------------|--|
| NBP2-35179-5ug | Recombinant Mouse CCL21/6Ckine Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| 457-6C-025 | CCL21/6Ckine [Unconjugated] |
| 485-MI-100 | IFN-gamma [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB457I

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

