

Product Datasheet

CD47 Antibody (472603) [FITC] FAB4670F

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB4670F

Updated 8/9/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB4670F



FAB4670F

CD47 Antibody (472603) [FITC]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	472603
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	FITC
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	961
Gene Symbol	CD47
Species	Human
Specificity/Sensitivity	Detects human CD47. Stains human CD47 transfectants but not irrelevant transfectants.
Immunogen	Mouse myeloma cell line NS0-derived recombinant human CD47 Gln19-Pro139 Accession # Q08722

Product Application Details	
Applications	Flow Cytometry, Blockade of Receptor-ligand Interaction, CyTOF-ready
Recommended Dilutions	Flow Cytometry, Blockade of Receptor-ligand Interaction, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

CD47 Antibody (472603) [FITC] [FAB4670F] - Vial of FITC conjugated antibody. FITC is optimally excited at 498 nm by the Blue laser (488 nm) and has an emission maximum of 519 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB4670F

NBP1-96849	Mouse IgG1 Isotype Control (MG1) [FITC]
NBP2-52320-0.05mg	Recombinant Human CD47 hIgG-His Protein
210-TA-005	TNF-alpha [Unconjugated]
1866-CD-050	CD47 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB4670F

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

