

# Product Datasheet

## PKC epsilon Antibody (666823) [FITC] FAB5134F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB5134F](http://www.novusbio.com/FAB5134F)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB5134F](http://www.novusbio.com/reviews/destination/FAB5134F)



**FAB5134F**

PKC epsilon Antibody (666823) [FITC]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	666823
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	FITC
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	PBS
Product Description	
Host	Mouse
Gene ID	5581
Gene Symbol	PRKCE
Species	Human
Specificity/Sensitivity	Detects human PKC $\epsilon$ in direct ELISAs. In direct ELISAs, less than 20% cross-reactivity with recombinant human (rh) PKC $\delta$ and rhPKC $\iota$ is observed. No cross-reactivity with rhPKC $\alpha$ , $\beta$ 1, $\beta$ 2, $\gamma$ , $\eta$ , $\mu$ , $\nu$ , $\theta$ , or $\zeta$ is observed.
Immunogen	<i>E. coli</i> -derived recombinant human PKC $\epsilon$ Gln580-Pro737 Accession # Q02156
Product Application Details	
Applications	Immunohistochemistry
Recommended Dilutions	Immunohistochemistry
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB5134F**

---

NBP2-27229	Mouse IgG2b Isotype Control (MPC-11) [FITC]
NBP2-54970PEP	PKC epsilon Recombinant Protein Antigen
236-EG-200	EGF [Unconjugated]
NBL1-14774	PKC epsilon Overexpression Lysate

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB5134F](http://www.novusbio.com/reviews/submit/FAB5134F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

