

# Product Datasheet

## ROR gamma t/RORC2/NR1F3 Antibody (1181A) [PE] FAB61093P

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB61093P](http://www.novusbio.com/FAB61093P)

Updated 10/15/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB61093P](http://www.novusbio.com/reviews/destination/FAB61093P)



**FAB61093P**

ROR gamma t/RORC2/NR1F3 Antibody (1181A) [PE]


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	1181A
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	PE
Purity	Protein A or G purified from cell culture supernatant
Buffer	PBS

Product Description	
Host	Rabbit
Gene ID	6097
Species	Human
Specificity/Sensitivity	Detects recombinant human ROR $\gamma$ t/RORC2/NR1F3 by direct ELISAs.
Immunogen	<i>E.coli</i> derived recombinant human ROR $\gamma$ T/RORC2/NR1F3 Met1-Gln100 Accession # P51449

Product Application Details	
Applications	Intracellular Staining by Flow Cytometry
Recommended Dilutions	Intracellular Staining by Flow Cytometry
Application Notes	Optimal dilution of this antibody should be experimentally determined.

**Images**

ROR gamma t/RORC2/NR1F3 Antibody (1181A) [PE] [FAB61093P] - Vial of PE conjugated antibody. PE has two excitation maxima, 498 nm excited by the Blue laser (488 nm) and 565 nm excited by the Yellow-Green laser (561 nm). Both result in emission at 578 nm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB61093P**

---

NBP2-24983	Rabbit IgG Isotype Control [PE]
M6000B-1	IL-6 [HRP]
DY421-05	IL-17/IL-17A [Biotin]
NB100-39002	FoxP3 Antibody - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB61093P](http://www.novusbio.com/reviews/submit/FAB61093P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

