

# Product Datasheet

## CD23/Fc epsilon RII Antibody (691632) [HRP] FAB6900H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6900H](http://www.novusbio.com/FAB6900H)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6900H](http://www.novusbio.com/reviews/destination/FAB6900H)



**FAB6900H**

CD23/Fc epsilon RII Antibody (691632) [HRP]

| Product Information         |                                                                                                                                                  |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                           |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                      |
| Storage                     | Store at 4C in the dark.                                                                                                                         |
| Clonality                   | Monoclonal                                                                                                                                       |
| Clone                       | 691632                                                                                                                                           |
| Preservative                | No Preservative                                                                                                                                  |
| Isotype                     | IgG2a                                                                                                                                            |
| Conjugate                   | HRP                                                                                                                                              |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                       |
| Buffer                      | PBS                                                                                                                                              |
| Product Description         |                                                                                                                                                  |
| Host                        | Rat                                                                                                                                              |
| Gene ID                     | 2208                                                                                                                                             |
| Gene Symbol                 | Fcer2a                                                                                                                                           |
| Species                     | Mouse                                                                                                                                            |
| Specificity/Sensitivity     | Detects mouse CD23/Fcε RII in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human (rh) Fcε RIA or rhFcε RIG is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant mouse CD23/Fcε RII Glu50-Pro331<br>Accession # P20693                                            |
| Product Application Details |                                                                                                                                                  |
| Applications                | Flow Cytometry, CyTOF-ready                                                                                                                      |
| Recommended Dilutions       | Flow Cytometry, CyTOF-ready                                                                                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                           |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB6900H**

---

|                |                                               |
|----------------|-----------------------------------------------|
| NBP2-35242-5ug | Recombinant Human CD23/Fc epsilon RII Protein |
| 210-TA-005     | TNF-alpha [Unconjugated]                      |
| DCD230         | CD23/Fc epsilon RII [HRP]                     |
| M6000B-1       | IL-6 [HRP]                                    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6900H](http://www.novusbio.com/reviews/submit/FAB6900H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

