

Product Datasheet

4-1BB/TNFRSF9/CD137 Antibody (2356C) [PE/Cy5.5] FAB8382PECY55

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB8382PECY55

Updated 2/24/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB8382PECY55



FAB8382PECY55

4-1BB/TNFRSF9/CD137 Antibody (2356C) [PE/Cy5.5]

| Product Information | |
|----------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. Do not freeze. |
| Clonality | Monoclonal |
| Clone | 2356C |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | PE/Cy5.5 |
| Purity | Protein A or G purified from cell culture supernatant |
| Buffer | PBS |

| Product Description | |
|--------------------------------|--|
| Host | Rabbit |
| Gene ID | 3604 |
| Species | Human |
| Specificity/Sensitivity | Detects human 4-1BB/TNFRSF9/CD137 in direct ELISAs. |
| Immunogen | Chinese Hamster Ovary cell line CHO-derived human 4-1BB/TNFRSF9/CD137 Leu24-His183 Accession # Q07011 |

| Product Application Details | |
|------------------------------------|--|
| Applications | Flow Cytometry |
| Recommended Dilutions | Flow Cytometry |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. For optimal results using our Tandem dyes, please avoid prolonged exposure to light or extreme temperature fluctuations. These can lead to irreversible degradation or decoupling. When staining intracellular targets, specific attention to the fixation and permeabilization steps in your flow protocol may be required. Please contact our technical support team at technical@novusbio.com if you have any questions. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB8382PECY55

| | |
|------------------|---|
| NBP2-24891PECY55 | Rabbit IgG Isotype Control [PE/Cy5.5] |
| NBP2-34958-5ug | Recombinant Human 4-1BB/TNFRSF9/CD137 Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| 937-4B-050 | 4-1BB/TNFRSF9/CD137 |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB8382PECY55

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

