

# Product Datasheet

## **BID Antibody (91508) [DyLight 350] FAB860D**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB860D](http://www.novusbio.com/FAB860D)

Updated 5/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB860D](http://www.novusbio.com/reviews/destination/FAB860D)



**FAB860D**

BID Antibody (91508) [DyLight 350]

| Product Information         |                                                                                                                                                                                                                                                                                                                     |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                              |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                         |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                            |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                          |
| Clone                       | 91508                                                                                                                                                                                                                                                                                                               |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                  |
| Isotype                     | IgG1                                                                                                                                                                                                                                                                                                                |
| Conjugate                   | DyLight 350                                                                                                                                                                                                                                                                                                         |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                          |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                  |
| Product Description         |                                                                                                                                                                                                                                                                                                                     |
| Host                        | Rat                                                                                                                                                                                                                                                                                                                 |
| Gene ID                     | 637                                                                                                                                                                                                                                                                                                                 |
| Gene Symbol                 | BID                                                                                                                                                                                                                                                                                                                 |
| Species                     | Human, Mouse                                                                                                                                                                                                                                                                                                        |
| Specificity/Sensitivity     | Detects human and mouse BID in Western blots. Detects full-length mouse BID and the 15 kDa carboxyl-terminal fragment (tBID) generated by cleavage with Caspase-8. The antibody also detects human BID. However, in Western blots of recombinant human and mouse BID, the antibody is most sensitive for mouse BID. |
| Immunogen                   | <i>E. coli</i> -derived recombinant mouse BID<br>Met1-Asp195<br>Accession # P70444                                                                                                                                                                                                                                  |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                   |
| Product Application Details |                                                                                                                                                                                                                                                                                                                     |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                        |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                                        |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                              |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB860D**

---

|               |                                            |
|---------------|--------------------------------------------|
| NBP1-86187PEP | BID Recombinant Protein Antigen            |
| 210-TA-005    | TNF-alpha [Unconjugated]                   |
| 882-B8-050    | BID [Unconjugated]                         |
| AF835         | Caspase-3 Antibody [Unconjugated] - Active |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB860D](http://www.novusbio.com/reviews/submit/FAB860D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

