

Product Datasheet

Smad2/3 [p Ser467, p Ser465] Antibody (1074A) [DyLight 405]

FAB8935E

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB8935E

Updated 5/13/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB8935E



FAB8935E

Smad2/3 [p Ser467, p Ser465] Antibody (1074A) [DyLight 405]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	1074A
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 405
Purity	Protein A or G purified from cell culture supernatant
Buffer	50mM Sodium Borate

Product Description	
Host	Rabbit
Gene ID	4087
Species	Human
Specificity/Sensitivity	Detects human Smad2 and Smad3 when dually phosphorylated at S465 and S467 in Western blots.
Immunogen	Phosphopeptide containing the human Smad2 S465/S467 site Accession # NP_005892
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB8935E

NBP2-24891V	Rabbit IgG Isotype Control [DyLight 405]
7754-BH-005/CF	TGF-beta 1 [Unconjugated]
AF467	EphB2 Antibody [Unconjugated]
MAB2029	Smad7 Antibody (293039) [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB8935E

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

