

# Product Datasheet

## **p70 S6 Kinase/S6K [p Ser424, p Thr421] Antibody (1143C) [PerCP] FAB8965C**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB8965C](http://www.novusbio.com/FAB8965C)

Updated 5/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB8965C](http://www.novusbio.com/reviews/destination/FAB8965C)



**FAB8965C**

p70 S6 Kinase/S6K [p Ser424, p Thr421] Antibody (1143C) [PerCP]

| <b>Product Information</b> |                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 0.1 ml                                                                                      |
| <b>Concentration</b>       | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>             | Store at 4C in the dark.                                                                    |
| <b>Clonality</b>           | Monoclonal                                                                                  |
| <b>Clone</b>               | 1143C                                                                                       |
| <b>Preservative</b>        | 0.05% Sodium Azide                                                                          |
| <b>Isotype</b>             | IgG                                                                                         |
| <b>Conjugate</b>           | PerCP                                                                                       |
| <b>Purity</b>              | Protein A or G purified from cell culture supernatant                                       |
| <b>Buffer</b>              | PBS                                                                                         |

| <b>Product Description</b>     |                                                                  |
|--------------------------------|------------------------------------------------------------------|
| <b>Host</b>                    | Rabbit                                                           |
| <b>Gene ID</b>                 | 6198                                                             |
| <b>Gene Symbol</b>             | RPS6KB1                                                          |
| <b>Species</b>                 | Human, Mouse, Rat                                                |
| <b>Specificity/Sensitivity</b> | Detects human p70 S6 Kinase when phosphorylated at T421/S424.    |
| <b>Immunogen</b>               | Phosphopeptide containing the human p70 S6 Kinase T421/S424 site |

| <b>Product Application Details</b> |                                                                        |
|------------------------------------|------------------------------------------------------------------------|
| <b>Applications</b>                | Western Blot                                                           |
| <b>Recommended Dilutions</b>       | Western Blot                                                           |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB8965C**

---

|               |                                             |
|---------------|---------------------------------------------|
| NBP2-24891PCP | Rabbit IgG Isotype Control [PerCP]          |
| NBP2-23649PEP | p70 S6 Kinase/S6K Antibody Blocking Peptide |
| 291-G1-200    | IGF-I/IGF-1 [Unconjugated]                  |
| 896-KS-010    | p70 S6 Kinase/S6K [Unconjugated]            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB8965C](http://www.novusbio.com/reviews/submit/FAB8965C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

