

# Product Datasheet

## EpCAM/TROP1 Antibody (G8.8R) [Janelia Fluor® 669] FAB8998JF669

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB8998JF669](http://www.novusbio.com/FAB8998JF669)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB8998JF669](http://www.novusbio.com/reviews/destination/FAB8998JF669)



**FAB8998JF669**

EpCAM/TROP1 Antibody (G8.8R) [Janelia Fluor® 669]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | G8.8R                                                                                       |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG2a                                                                                       |
| Conjugate           | Janelia Fluor 669                                                                           |
| Purity              | 0                                                                                           |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                       |
|-------------------------|---------------------------------------------------------------------------------------|
| Host                    | Rat                                                                                   |
| Gene ID                 | 4072                                                                                  |
| Gene Symbol             | EPCAM                                                                                 |
| Species                 | Mouse                                                                                 |
| Specificity/Sensitivity | Detects mouse EpCAM/TROP-1.                                                           |
| Immunogen               | TE-71 Thymic epithelial cell line                                                     |
| Notes                   | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus. |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Flow Cytometry, CyTOF-reported                                         |
| Recommended Dilutions       | Flow Cytometry, CyTOF-reported                                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB8998JF669**

---

|                   |                                           |
|-------------------|-------------------------------------------|
| NBP2-52190-0.05mg | Recombinant Human EpCAM/TROP1 His Protein |
| 202-IL-010        | IL-2 [Unconjugated]                       |
| 960-EP-050        | EpCAM/TROP1 [Unconjugated]                |
| AF748             | E-Cadherin Antibody [Unconjugated]        |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB8998JF669](http://www.novusbio.com/reviews/submit/FAB8998JF669)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

