

# Product Datasheet

## IL-10 Antibody (2050B) [Allophycocyanin] FAB9210A

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB9210A](http://www.novusbio.com/FAB9210A)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB9210A](http://www.novusbio.com/reviews/destination/FAB9210A)



**FAB9210A**

IL-10 Antibody (2050B) [Allophycocyanin]

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C in the dark.                                                                    |
| <b>Clonality</b>     | Monoclonal                                                                                  |
| <b>Clone</b>         | 2050B                                                                                       |
| <b>Preservative</b>  | 0.05% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG2b                                                                                       |
| <b>Conjugate</b>     | Allophycocyanin                                                                             |
| <b>Purity</b>        | Protein A or G purified from cell culture supernatant                                       |
| <b>Buffer</b>        | PBS                                                                                         |

**Product Description**

|                                |                                                                                                                       |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>                    | Rabbit                                                                                                                |
| <b>Gene ID</b>                 | 3586                                                                                                                  |
| <b>Species</b>                 | Human                                                                                                                 |
| <b>Specificity/Sensitivity</b> | Detects human IL-10 in direct ELISAs.                                                                                 |
| <b>Immunogen</b>               | <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant human IL-10 Met1-Asn178<br>Accession # P22301 |

**Product Application Details**

|                              |                                                                        |
|------------------------------|------------------------------------------------------------------------|
| <b>Applications</b>          | Intracellular Staining by Flow Cytometry, Neutralization               |
| <b>Recommended Dilutions</b> | Neutralization, Intracellular Staining by Flow Cytometry               |
| <b>Application Notes</b>     | Optimal dilution of this antibody should be experimentally determined. |

**Images**

IL-10 Antibody (2050B) [Allophycocyanin] [FAB9210A] - Vial of APC conjugated antibody. APC is optimally excited at 650 nm by the Red laser (633 or 640 nm) and has an emission maximum of 660 nm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB9210A**

---

|                 |                                 |
|-----------------|---------------------------------|
| NBP2-35101-10ug | Recombinant Mouse IL-10 Protein |
| 210-TA-005      | TNF-alpha [Unconjugated]        |
| DY417-05        | IL-10 [Biotin]                  |
| M6000B-1        | IL-6 [HRP]                      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB9210A](http://www.novusbio.com/reviews/submit/FAB9210A)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

