

# Product Datasheet

## CD30/TNFRSF8 Antibody (Hu8) [PE] FAB9576P

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB9576P](http://www.novusbio.com/FAB9576P)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB9576P](http://www.novusbio.com/reviews/destination/FAB9576P)



**FAB9576P****CD30/TNFRSF8 Antibody (Hu8) [PE]**

| <b>Product Information</b>         |                                                                                                                                                                                                   |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                            |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                       |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                          |
| <b>Clonality</b>                   | Monoclonal                                                                                                                                                                                        |
| <b>Clone</b>                       | Hu8                                                                                                                                                                                               |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                |
| <b>Isotype</b>                     | IgG1                                                                                                                                                                                              |
| <b>Conjugate</b>                   | PE                                                                                                                                                                                                |
| <b>Purity</b>                      | Protein A or G purified from cell culture supernatant                                                                                                                                             |
| <b>Buffer</b>                      | PBS                                                                                                                                                                                               |
| <b>Product Description</b>         |                                                                                                                                                                                                   |
| <b>Host</b>                        | Human                                                                                                                                                                                             |
| <b>Gene ID</b>                     | 943                                                                                                                                                                                               |
| <b>Species</b>                     | Human                                                                                                                                                                                             |
| <b>Specificity/Sensitivity</b>     | Detects human CD30/TNFRSF8 in direct ELISAs. This non-therapeutic antibody uses the same variable region sequence as the therapeutic antibody Brentuximab. This product is for research use only. |
| <b>Immunogen</b>                   | Human CD30/TNFRSF8                                                                                                                                                                                |
| <b>Product Application Details</b> |                                                                                                                                                                                                   |
| <b>Applications</b>                | Flow Cytometry, CyTOF-ready                                                                                                                                                                       |
| <b>Recommended Dilutions</b>       | Flow Cytometry, CyTOF-ready                                                                                                                                                                       |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB9576P**

---

|                |                                            |
|----------------|--------------------------------------------|
| NBP2-22660     | Recombinant Human CD30/TNFRSF8 His Protein |
| 210-TA-005     | TNF-alpha [Unconjugated]                   |
| 813-CD-100     | CD30/TNFRSF8 [Unconjugated]                |
| 6507-IL-010/CF | IL-4 [Unconjugated]                        |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB9576P](http://www.novusbio.com/reviews/submit/FAB9576P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

