

Product Datasheet

MVP Antibody (2H3-1A6) - Azide and BSA Free H00009961-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00009961-M01

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00009961-M01



H00009961-M01

MVP Antibody (2H3-1A6) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2H3-1A6
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4

Product Description	
Description	Novus Biologicals Mouse MVP Antibody (2H3-1A6) - Azide and BSA Free (H00009961-M01) is a monoclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. Anti-MVP Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	9961
Gene Symbol	MVP
Species	Human
Specificity/Sensitivity	MVP - major vault protein
Immunogen	MVP (AAH15623, 1 a.a. ~ 893 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MATEEFIIRIPPYHYIHVLDQNSNVS RVEVGPKTYIRQDNERVLFAPMRMVTVPP RHYCTVANPVS RDAQGLV LFDVTGQVRLRHADLEIRLAQDPFPLYPGEVLEKDI TPLQVVLPNTALHLKALLDFEDKDGKVVAGDEWLFEGPGTYIPRKEVEVVEII QATIIRQNQALRLRARKECWDRD GKERVTGEEWLVTTVGAYLPAVFEEVLDLV DAVILTEKTALHLRARRNFRDFR GVSRRTGEEWLVTVQDTEAHVPDVHEEVLG VVPITTLGPHNYCVILDPVGP DGNKQLGQKRVVKGEKSFFLQPGEQLEQGIQD VYVLSEQQGLLLRALQPLEEGEDEEKVSHQAGDHWLIRGPLEYVPSAKVEVVE ERQAIPLDENEGIYVQDVKTGKVRVAVIGSTYMLTQDEVLWEKELPPGVEELLNK GQDPLADRGEKDTAKSLQPLAPRNKTRVVS YRVPHNAAVQVYDYREKRARVV FGPELVSLGPEEQFTVLSLSAGRPKRPHARRALCLLLGPDFFTDVITETADHAR LQLQLAYNWHFEVNDRKDPQETAKLFSVPDFVGDACKAIASRVRGAVASVTFD DFHKNSARIIRTAVFGFETSEAKGPDGMALPRPRDQAVFPQNGLVSSVDVQS VEPVDQRTRDALQRSVQLAIEITTNSQEAAAKHEAQRLEQEARGRLERQKILDQ SEAEKARKELLELEALSMAVESTGTAKAEAESRAEAARIEGEGSVLQAKLKAQA LAIETEAE LQRVQKVRELELVYARAQLELEVSKAQQLAEVEVKKFKQMTEAIGP STIRDLAVAGPEMQVKLLQSLGLKSTLITDGSTPINLFNTAFGLLGMGPEGQPL GRRVASGSPGEGISPQSAQAPQAPGDNHVVPVLR
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry




Recommended Dilutions	Western Blot 1:500, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin
Application Notes	Antibody reactive against cell lysate and recombinant protein for Western Blot. Has also been used for immunofluorescence, immunohistochemistry (paraffin), and ELISA.

Images


Western Blot: MVP Antibody (2H3-1A6) [H00009961-M01] - Western blot of several proteins that were increased or decreased in DBA patients including dysferlin (DYSF), gamma-adducin (ADD3), major vault protein (MVP), serine palmitoyltransferase 1 and 2 (SPTLC1 and SPTLC2), HLA-ABC, and intercellular adhesion molecule 4 (ICAM4) with b-actin (ACTB) used as a loading control. Six individual control samples that were included in the pooled controls for proteomic analysis were used for comparison to individual DBA patient samples. Image collected and cropped by CiteAb from the following publication (dx.plos.org/10.1371/journal.pone.0085504), licensed under a CC-BY license.




Immunocytochemistry/Immunofluorescence: MVP Antibody (2H3-1A6) [H00009961-M01] - Analysis of monoclonal antibody to MVP on HeLa cell. Antibody concentration 10 ug/ml.




Immunohistochemistry-Paraffin: MVP Antibody (2H3-1A6) [H00009961-M01] - Analysis of monoclonal antibody to MVP on formalin-fixed paraffin-embedded human lung, adenocarcinoma. Antibody concentration 3 ug/ml.




Western Blot: MVP Antibody (2H3-1A6) [H00009961-M01] - MVP monoclonal antibody (M01), clone 2H3-1A6 Analysis of MVP expression in Hela S3 NE.



ELISA: MVP Antibody (2H3-1A6) [H00009961-M01] - Detection limit for recombinant GST tagged MVP is approximately 1ng/ml as a capture antibody.



Western Blot: MVP Antibody (2H3-1A6) [H00009961-M01] - Immunoblot validation & characterization of dysferlin. A) Western blot of several proteins that were increased or decreased in DBA patients including dysferlin (DYSF), gamma-adducin (ADD3), major vault protein (MVP), serine palmitoyltransferase 1 & 2 (SPTLC1 & SPTLC2), HLA-ABC, & intercellular adhesion molecule 4 (ICAM4) with β -actin (ACTB) used as a loading control. Six individual control samples that were included in the pooled controls for proteomic analysis were used for comparison to individual DBA patient samples. B) Confocal images of intact RBCs from DBA patient D1 & healthy control C2 showing localization of dysferlin (green) at the membrane where spectrin (red) is localized. Dysferlin fluorescence signal in the membrane of DBA RBCs is present in punctuate patterns & is evident above the autofluorescent background seen in the cytoplasm as highlighted in the zoomed images on the right panel. C) Western blot showing the enrichment of dysferlin in the RBC ghosts (WG) & the inside-out vesicles (IOV) but not in the RBC membrane cytoskeletal extract (CE) of DBA patients D1 & D2. D) Distribution of dysferlin spectral counts versus gel slice number for DBA samples D1 & D2. Higher abundance of dysferlin in D1 with high molecular weight dysferlin crosslinks & low mass protein fragments correlate with immunoblots shown in panel A. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/24454878>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Rahimi AS, Wilson DD, Saylor DK et al. p16, Cyclin D1, and HIF-1alpha Predict Outcomes of Patients with Oropharyngeal Squamous Cell Carcinoma Treated with Definitive Intensity-Modulated Radiation Therapy. *Int J Otolaryngol.* 2012-01-01 [PMID: 22888357] (IF/IHC)

Pesciotta EN, Sriswasdi S, Tang HY et al. Dysferlin and other non-red cell proteins accumulate in the red cell membrane of Diamond-Blackfan Anemia patients. *PLoS ONE* 2014-01-23 [PMID: 24454878] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00009961-M01

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00009961-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

