

Product Datasheet

FGFR1OP Antibody (2B1) - Azide and BSA Free H00011116-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 19


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00011116-M01

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00011116-M01




H00011116-M01**FGFR1OP Antibody (2B1) - Azide and BSA Free**

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2B1
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	IgG purified
Buffer	PBS (pH 7.4)
Product Description	
Description	Novus Biologicals Mouse FGFR1OP Antibody (2B1) - Azide and BSA Free (H00011116-M01) is a monoclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. Anti-FGFR1OP Antibody: Cited in 19 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	11116
Gene Symbol	CEP43
Species	Human, Mouse, Rabbit
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 28625565). Other species not tested. Rabbit reactivity reported from a verified customer review.
Specificity/Sensitivity	FGFR1OP - FGFR1 oncogene partner
Immunogen	FGFR1OP (AAH11902, 1 a.a ~ 379 a.a) full length recombinant protein with GST tag.MAATAAAVVAEEDTELRLDQVTLNSGVLNRIKAE LRAAVFLALEEQEKVE NKTPLVNESLRKFLNTKDGRLVASLVAEFLQFFNLDFTLAVFQPETSTLQGLEGR ENLARDLGIIIEAGTVGGPLLLEVIRRCQQKEKGPTTGEGALDLSDVHSPPKS PEGKTSAQTTSPSKKANDEANQSDTSVSLSEPKSKSSLHLLSHETKIGSFSLNRT LDGKDKAGLCPDEDDMEGDSFFDDPIPKPEKTYGLRNEPRK
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Western Blot 1:500, ELISA 1:100 - 1:2000, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:10 - 1:500, Immunohistochemistry-Frozen
Application Notes	Antibody reactive against recombinant protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control. It is also useful for immunofluorescence. This FGFR1OP Antibody (2B1) is validated for IHC-Fr from a verified customer review.



Images


Western Blot: FGFR1OP Antibody (2B1) [H00011116-M01] - Analysis of FGFR1OP expression in transfected 293T cell line by FGFR1OP monoclonal antibody (M01), clone 2B1. Lane 1: FGFR1OP transfected lysate(43.1 KDa). Lane 2: Non-transfected lysate.




Immunocytochemistry/Immunofluorescence: FGFR1OP Antibody (2B1) [H00011116-M01] - Analysis of monoclonal antibody to FGFR1OP on HeLa cell. Antibody concentration 10 ug/mL.



Immunohistochemistry-Frozen: FGFR1OP Antibody (2B1) [H00011116-M01] - Rabbit frozen 4% paraformalehyde fixed brain sections. Double labeling with FGFR1OP (Basal bodies in green) and GT335 (cilia in red) antibodies of ependymal whole mount from brain of rabbit kit. IHC-Fr image submitted by a verified customer review.



ELISA: FGFR1OP Antibody (2B1) [H00011116-M01] - Detection limit for recombinant GST tagged FGFR1OP is approximately 0.03 ng/mL as a capture antibody.



Publications

- Wu CT, Hilgendorf KI, Bevacqua RJ et al. Discovery of ciliary G protein-coupled receptors regulating pancreatic islet insulin and glucagon secretion *Genes & Development* 2021-09-01 [PMID: 34385262]
- Gonçalves AB, Hasselbalch SK, Joensen BB et al. CEP78 functions downstream of CEP350 to control biogenesis of primary cilia by negatively regulating CP110 levels *eLife* 2021-07-14 [PMID: 34259627]
- Basso M, Mahuzier A, Ali SK, Marty A et al. Actin-based deformations of the nucleus control mouse multiciliated ependymal cell differentiation *Dev Cell* 2024-12-11 [PMID: 39662468]
- Mazo G, Soplop N, Wang WJ et al. Spatial Control of Primary Ciliogenesis by Subdistal Appendages Alters Sensation-Associated Properties of Cilia *Developmental Cell* 2016-11-01 [PMID: 27818179]
- Adelaida P, Keren H, Ann Y et al. Primary cilia on muscle stem cells are critical to maintain regenerative capacity and are lost during aging. *Nat Commun.* 2022-03-17 [PMID: 35301320]
- Racha C, Adham S, Xavier P et al. A Dual Protein-mRNA Localization Screen Reveals Compartmentalized Translation and Widespread Co-translational RNA Targeting. *Dev Cell.* 2020-08-04 [PMID: 32783880]
- Sayaka F, Hantian Q, Shohei N et al. ARL3 and ARL13B GTPases participate in distinct steps of INPP5E targeting to the ciliary membrane. *Biol Open.* 2021-09-28 [PMID: 34447983]
- Lena B, Anca M, Camille B et al. LRP2 contributes to planar cell polarity-dependent coordination of motile cilia function. *Cell Tissue Res.* 2023-02-11 [PMID: 36764939]
- Martina B, Nicolas K, Serena B et al. mTOR and S6K1 Drive Polycystic Kidney by the Control of Afadin-dependent Oriented Cell Division. *Nat Commun.* 2020-01-24 [PMID: 32581239]
- Tomoharu K, Roy N, Keene A et al. Myristoylated Neuronal Calcium Sensor-1 captures the ciliary vesicle at distal appendages. *bioRxiv.* 2023-01-10 [PMID: 36712037]
- Tomoharu K, Julia L, Saxton F et al. A hierarchical pathway for assembly of the distal appendages that organize primary cilia. *bioRxiv.* 2023-01-10 [PMID: 36711481]
- Takuya K, Yamato I, Tomoaki H et al. Cooperation of the IFT-A complex with the IFT-B complex is required for ciliary retrograde protein trafficking and GPCR import. *Mol Biol Cell.* 2020-11-11 [PMID: 33175651]
- More publications at <http://www.novusbio.com/H00011116-M01>





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00011116-M01

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00011116-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

