

Product Datasheet

AIPL1 Antibody (3A3) - Azide and BSA Free H00023746-M23

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00023746-M23

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00023746-M23



H00023746-M23


AIPL1 Antibody (3A3) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3A3
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Novus Biologicals Mouse AIPL1 Antibody (3A3) - Azide and BSA Free (H00023746-M23) is a monoclonal antibody validated for use in ELISA and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	23746
Gene Symbol	AIPL1
Species	Human
Specificity/Sensitivity	AIPL1 - aryl hydrocarbon receptor interacting protein-like 1 (3A3)
Immunogen	AIPL1 (NP_055151.3, 1 a.a. ~ 101 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MDAALLLNVEGVKKTILHGGTGELPNFITGSRVIFHFRTMKCDEERTVIDDSRQV GQPMHIIIIGNMFKLEVWEILLTSMRVHEVAEFWCDTIHTGVYPILS
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	ELISA, Immunoprecipitation, Sandwich ELISA
Recommended Dilutions	ELISA, Immunoprecipitation, Sandwich ELISA
Application Notes	This product is useful for ELISA.



Images

Sandwich ELISA: AIPL1 Antibody (3A3) [H00023746-M23] - Detection limit for recombinant GST tagged AIPL1 is 0.3 ng/ml as a capture antibody.



Immunoprecipitation: AIPL1 Antibody (3A3) [H00023746-M23] - Analysis of AIPL1 transfected lysate using anti-AIPL1 monoclonal antibody and Protein A Magnetic Bead, and immunoblotted with AIPL1 rabbit polyclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00023746-M23

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
NBP2-55506PEP	AIPL1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00023746-M23

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

