

Product Datasheet

SAMM50 Antibody (2A9) - Azide and BSA Free H00025813-M04

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00025813-M04

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00025813-M04



H00025813-M04**SAMM50 Antibody (2A9) - Azide and BSA Free**


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2A9
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4

Product Description	
Description	Novus Biologicals Mouse SAMM50 Antibody (2A9) - Azide and BSA Free (H00025813-M04) is a monoclonal antibody validated for use in WB and ELISA. Anti-SAMM50 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	25813
Gene Symbol	SAMM50
Species	Human, Rat
Reactivity Notes	Other species not tested.
Specificity/Sensitivity	SAMM50 (2A9)
Immunogen	SAMM50 (AAH07830, 1 a.a. ~ 468 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MGTVHARSLEPLPSSGPDFGGLGEEAEFVEVEPEAKQEILENKDVVVQHVHFD GLGRTKDDIIICEIGDVFKAKNLIEVMRKSHEAREKLLRLGIFRQVDVLIDTCQGD DALPNGLDVTFEVTELRRLTGSYNTMVGNNEGSMVLGLKLPNLLGRAEKVTFQ FSYGTKETSYGLSFFKPRPGNFERNFVSVNLYKVTGQFPWSSLRETDGRGMSAE YSFPIWKTSHTVKWEVWRELGCLSRASFAVRKESGHSLKSSLSHAMVIDSR NSSILPRRGALLKVNQELAGYTGGDVSEFIKEDFELQLNKQLIFDSVFSASFVGG MLVPIGDKPSSIADRFYLGGPSTVRGFSMHSIGPQSEG DYLGGEAYWAGLHLY TPLPFRPGQGGFGE LFRTHFFLNAGNLCNLNYGEGPKAHIRKLAECIRWSYGA GIVLRLGNIARLELNVCVPMGVQTGDRICDGVQFGAGIRFL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500, ELISA 1:100-1:2000
Application Notes	Antibody reactive against cell lysate and recombinant protein for western blot. It has also been used for ELISA.

Images

Western Blot: SAMM50 Antibody (2A9) [H00025813-M04] - Analysis of SAMM50 expression in A-431 (Cat # L015V1).



Publications

Joseph AM, Hood DA. Plasticity of TOM complex assembly in skeletal muscle mitochondria in response to chronic contractile activity. *Mitochondrion*. 2011-11-25 [PMID: 22142511]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00025813-M04

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00025813-M04

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

