

Product Datasheet

PI4K2A Antibody - Azide and BSA Free H00055361-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00055361-B01P

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00055361-B01P



H00055361-B01P

PI4K2A Antibody - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	IgG purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Mouse PI4K2A Antibody - Azide and BSA Free (H00055361-B01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	55361
Gene Symbol	PI4K2A
Species	Human
Specificity/Sensitivity	PI4KII - phosphatidylinositol 4-kinase type II,
Immunogen	PI4KII (NP_060895.1, 1 a.a. - 479 a.a.) full-length human protein. MDETSPLVSPERAQPPDYTFPSGSGAHFPQVPGGAVRVAAAAGSGPSPPGSP GHDREERQPLLD RARGAAAQGGQTQTVAQAQALAAQAAAAHAAQAHRENE FPEDPEFEAVVRQAELAIERCIFPERIYQGSSGSYFVKDPQGRIIAVFKPKNEEP YGHLPKWKWLQKLCCPCCFGRDCLVLNQGYLSEAGASLVDQKLELNIVPR TKVVYLASETFNYSIDRVKSRGKRLALEKVPKVGQRFNRIGLPPKVGSGFQLFV EGYKDADYWLRRFEAEPLPENTNRQLLLQFERLVLDYIIRNTDRGNDNWLKY DCPMDSSSSRD TDWVVVKEPVKVAIDNGLAFPLKHPDSWRAYPFYWAWLP QAKVPFSQEIKDLILPKISDPNFVKDLEEDLYELFKKDPGFDRGQFHKQIAVMRG QILNLTQALKDNKSPLHLVQMPPVIVETARSHQRSSSESYSYTSFQSRKPFSSW W
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:500
Application Notes	Antibody reactive against transfected lysate in western blot. May also be used for ELISA.




Images

Western Blot: PI4K2A Antibody [H00055361-B01P] - Analysis of PI4K2A expression in HepG2.



Western Blot: PI4K2A Antibody [H00055361-B01P] - Analysis of PI4KII expression in transfected 293T cell line by PI4KII polyclonal antibody. Lane 1: PI4KII transfected lysate(52.69 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00055361-B01P

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97019-5mg	Mouse IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00055361-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

