

Product Datasheet

VPS35 Antibody (2D3) - Azide and BSA Free H00055737-M02

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 9


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00055737-M02

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00055737-M02



H00055737-M02

VPS35 Antibody (2D3) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2D3
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Novus Biologicals Mouse VPS35 Antibody (2D3) - Azide and BSA Free (H00055737-M02) is a monoclonal antibody validated for use in WB, ELISA, ICC/IF and IP. Anti-VPS35 Antibody: Cited in 9 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	55737
Gene Symbol	VPS35
Species	Human, Mouse, Rat
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 25190809). Other species not tested.
Specificity/Sensitivity	VPS35 (2D3)
Immunogen	VPS35 (NP_060676, 697 a.a. ~ 796 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MECLKKALKIANQCMDPSLQVQLFIEILNRYIYFYEKENDAVTIQVLNQLIQKIRE DLPNLESSEETE QINKHFHNTLEHLRLRRESPESSEGP IYEG LIL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500, ELISA 1:100-1:2000, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Application Notes	Antibody Reactive against Recombinant Protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control. Use in Immunocytochemistry/immunofluorescence reported in scientific literature (PMID 28146471).



Images


Western Blot: VPS35 Antibody (2D3) [H00055737-M02] - VPS35 monoclonal antibody (M02), clone 2D3. Analysis of VPS35 expression in human colon.



Western Blot: VPS35 Antibody (2D3) [H00055737-M02] - VPS35 monoclonal antibody (M02), clone 2D3. Analysis of VPS35 expression in HeLa.



ELISA: VPS35 Antibody (2D3) [H00055737-M02] - Detection limit for recombinant GST tagged VPS35 is 0.03 ng/ml as a capture antibody.



Publications

Abhishek M, Rodolfo Z, Vincent S et al. A new functional role of mitochondria-anchored protein ligase in peroxisome morphology in mammalian cells. *J Cell Biochem.* 2021-07-28 [PMID: 34322908]

Chelsie K, Anusha K, Austen M et al. Endosomal traffic and glutamate synapse activity are increased in VPS35 D620N mutant knock-in mouse neurons, and resistant to LRRK2 kinase inhibition. *Mol Brain.* 2021-09-16 [PMID: 34530877]

Follett J, Fox JD, Gustavsson EK et al. DNAJC13 p.Asn855Ser, implicated in familial parkinsonism, alters membrane dynamics of Sorting Nexin 1. *Neurosci Lett.* 2019-05-11 [PMID: 31082451]

Cataldi S, Follett J, Fox JD et al. Altered dopamine release and monoamine transporters in Vps35 p.D620N knock-in mice. *NPJ Parkinsons Dis* 2018-08-21 [PMID: 30155515]

Ayumu S, Sevan M, Prudent J, McBride H. Newly born peroxisomes are a hybrid of mitochondrial and ER-derived pre-peroxisomes *Nature.* [PMID: 28146471] (ICC/IF, Human)

Prosser DC, Tran D, Schooley A et al. A Novel, Retromer-Independent Role for Sorting Nexins 1 and 2 in RhoG-Dependent Membrane Remodeling. *Traffic.* 2010-08-03 [PMID: 20604901]

Munsie LN, Milnerwood AJ, Seibler P et al. Retromer-dependent neurotransmitter receptor trafficking to synapses is altered by the Parkinson's Disease VPS35 mutation p.D620N. *Hum Mol Genet.* 2014-11-21 [PMID: 25416282]

Brigidi GS, Santyr B, Shimell J et al. Activity-regulated trafficking of the palmitoyl-acyl transferase DHHC5. *Nat Commun.* 2015-09-03 [PMID: 26334723]

Kang Yunhee, Ge Yuan, Cassidy Robert M et al. A combined transgenic proteomic analysis and regulated trafficking of neuroligin-2. *J Biol Chem.* 2014-10-17 [PMID: 25190809] (Mouse)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00055737-M02

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00055737-M02

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

