

# Product Datasheet

## APMAP Antibody - Azide and BSA Free H00057136-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00057136-B01P](http://www.novusbio.com/H00057136-B01P)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00057136-B01P](http://www.novusbio.com/reviews/destination/H00057136-B01P)



**H00057136-B01P**

APMAP Antibody - Azide and BSA Free

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.05 mg                                                                                                          |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| Clonality           | Polyclonal                                                                                                       |
| Preservative        | No Preservative                                                                                                  |
| Isotype             | IgG                                                                                                              |
| Purity              | Protein A purified                                                                                               |
| Buffer              | PBS (pH 7.4)                                                                                                     |


| Product Description     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description             | Novus Biologicals Mouse APMAP Antibody - Azide and BSA Free (H00057136-B01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                                                                                                                                                                                                                                                   |
| Host                    | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Gene ID                 | 57136                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Gene Symbol             | APMAP                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Species                 | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Specificity/Sensitivity | C20orf3 - chromosome 20 open reading frame 3,                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Immunogen               | C20orf3 (NP_065392.1, 1 a.a. - 416 a.a.) full-length human protein.<br>MSEADGLRQRRPLRPQVVTDDDGGQAPEAKDGSSFSGRVFRVTFLMLAVSLTV<br>PLLGAMMLLESPIDPQPLSFKEPPLLLGVLHPNTKLRQAERLFENQLVGPESIAH<br>IGDVMFTGTADGRVVKLENGEIETIARFGSGPCKTRDDEPVCGRPLGIRAGPN<br>GTLFVADAYKGLFEVNPWKREVKLLLSSETPIEGKNMSFVNDLTVTQDGRKIYF<br>TDSSSKWQRRDYLLLVMEGTDDGRLLLEYDTVTREVKVLLDQLRFPNGVQLSP<br>AEDFVLVAETTMARIRRVYVSGLMKGGADLFVENMPGFPDNIRPSSSGGYWV<br>GMSTIRPNPGFSMLDFLSERPWIKRMIFKLFSQETVMKFVPRYSLVLELSDSGA<br>FRRSLHDPDGLVATYISEVHEHDGHLYLGSRSPFLCRLSLQAV |
| Notes                   | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

| Product Application Details |                                                                                                                                                                                   |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot                                                                                                                                                                      |
| Recommended Dilutions       | Western Blot 1:500                                                                                                                                                                |
| Application Notes           | Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control. |




## Images


Western Blot: APMAP Antibody [H00057136-B01P] - Analysis of C20orf3 expression in human liver.



Western Blot: APMAP Antibody [H00057136-B01P] - Analysis of C20orf3 expression in A-431.



Western Blot: APMAP Antibody [H00057136-B01P] - Analysis of C20orf3 expression in transfected 293T cell line by C20orf3 polyclonal antibody. Lane 1: C20orf3 transfected lysate(45.76 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to H00057136-B01P**

---

|                |                                                    |
|----------------|----------------------------------------------------|
| NBP2-33376H    | Blue Marker Antibody (6F4-F6) [HRP]                |
| HAF007         | Goat anti-Mouse IgG Secondary Antibody [HRP]       |
| NB7539         | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-97019-5mg | Mouse IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00057136-B01P](http://www.novusbio.com/reviews/submit/H00057136-B01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

