

Product Datasheet

Dynein intermediate chain 2 Antibody - Azide and BSA Free H00064446-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00064446-D01P

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00064446-D01P



H00064446-D01P**Dynein intermediate chain 2 Antibody - Azide and BSA Free**


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description	
Description	Novus Biologicals Rabbit Dynein intermediate chain 2 Antibody - Azide and BSA Free (H00064446-D01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	64446
Gene Symbol	DNAI2
Species	Human
Specificity/Sensitivity	DNAI2 - dynein, axonemal, intermediate chain 2,
Immunogen	DNAI2 (AAH39582.1, 1 a.a. - 593 a.a.) full-length human protein. MEIVYVYVKKRSEFGKQCNFSDRQAELNIDIMPNELAEQFVERNPVDTGIQCS ISMSEHEANSERFEMETRGVNHVEGGWPKDVNPLELEQTIRFRKKVEKDENVV NAIMQLGSIMEHCIKQNNADIYEEYFNDEEAMEVMEEDPSAKTINVFRDPQEIK RAATHLSWHPDGNRKLAVAYSCLDFQRAPVGMSSDSYIWDLENPNKPELALK PSSPLVTLEFNPKDSHVLLGGCYNGQIACWDTRKGLVAELSTIESSHRDPVY GTIWLQSKTGTECFASASTDGQVMWWDIRKMSEPTEVVILDITKKEQLENALGAI SLEFESTLPTKFMVGTEQGIVISCNRAKTSAEKIVCTFPGHHGPIYALQRNPFY PKNFLTVDWTARIWSEDSRESSIMWTKYHMAYLTDAAWSPVRPTVFFTTRM DGTLDIWDFMFEQCDPTLSLKDNGCLIACGSQGLGTTTLLEVSPGLSTLQRNEKN VASSMFERETRREKILEARHREMLRKEKGAEGRDEEQTDEELAVDLEALVSK AEEFFDIIFTELKKKEADAIKLTPVPPQPSPEEDQVVEEGEEAAGEEGDEEVE EDLA
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:500
Application Notes	Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control.

Images

Western Blot: Dynein intermediate chain 2 Antibody [H00064446-D01P] - Analysis of DNAI2 expression in transfected 293T cell line by DNAI2 polyclonal antibody. Lane 1: DNAI2 transfected lysate(67.50 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00064446-D01P

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00064446-D01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

