

Product Datasheet

Vasorin/SLIT-like 2 Antibody (4G7) - Azide and BSA Free H00114990-M05

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00114990-M05

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00114990-M05




H00114990-M05**Vasorin/SLIT-like 2 Antibody (4G7) - Azide and BSA Free**


Product Information	
Unit Size	0.05 mg
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4G7
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Ascites
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Novus Biologicals Mouse Vasorin/SLIT-like 2 Antibody (4G7) - Azide and BSA Free (H00114990-M05) is a monoclonal antibody validated for use in WB and ELISA. Anti-Vasorin/SLIT-like 2 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	114990
Gene Symbol	VASN
Species	Human
Specificity/Sensitivity	SLITL2 - slit-like 2 (Drosophila) (4G7)
Immunogen	VASN (NP_612449, 298 a.a. ~ 349 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. NPFNCVCPLSWFGPWVRESHVTLASPEETRCHFPKKNAGRLLLLLELDYADFGC
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500, ELISA
Application Notes	Antibody Reactive against Recombinant Protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control.

Images

Western Blot: Vasorin/SLIT-like 2 Antibody (4G7) [H00114990-M05] - Analysis of VASN expression in transfected 293T cell line by VASN monoclonal antibody (M05), clone 4G7. Lane 1: VASN transfected lysate (Predicted MW: 71.7 KDa). Lane 2: Non-transfected lysate.



ELISA: Vasorin/SLIT-like 2 Antibody (4G7) [H00114990-M05] - Detection limit for recombinant GST tagged SLITL2 is approximately 0.3ng/ml as a capture antibody.



Publications

Huang A, Dong J, Li S et al. Exosomal Transfer of Vasorin Expressed in Hepatocellular Carcinoma Cells Promotes Migration of Human Umbilical Vein Endothelial Cells. *Int. J. Biol. Sci.* 2015. 2015-01-01 [PMID: 26157350]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00114990-M05

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00114990-M05

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

