

Product Datasheet

GLIPR1L1 Antibody - Azide and BSA Free H00256710-B03P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00256710-B03P

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00256710-B03P



H00256710-B03P

GLIPR1L1 Antibody - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Mouse GLIPR1L1 Antibody - Azide and BSA Free (H00256710-B03P) is a polyclonal antibody validated for use in WB and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	256710
Gene Symbol	GLIPR1L1
Species	Human
Specificity/Sensitivity	Reacts with GLI pathogenesis-related 1 like 1.
Immunogen	GLIPR1L1 (NP_689992.1, 1 a.a. - 233 a.a.) full-length human protein. MALKNKFSCLWILGLCLVATTSSKIPSITDPHFIDNCIEAHNEWRGKVNPPAADM KYMIWDKGLAKMAKAWANQCKFEHNDCLDKSYKCYAAFEYVGENIWLGGIKS FTPRHAITAWYNETQFYDFDSLSCSRVCGHYTQLVWANSFYVGCAMCPNL GGASTAIFVCNYGPAGNFANMPPYVRGESCSSLCSKEEKCVKNLCKNPFLKPTG RAPQQTAFNPFSLGFLLRIF
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Flow Cytometry
Recommended Dilutions	Western Blot 1:500, Flow Cytometry
Application Notes	This antibody is reactive against transfected lysate in western blot, FACS analysis, and as a detection antibody in ELISA. It is also useful in FACS analysis.




Images

Western Blot: GLIPR1L1 Antibody [H00256710-B03P] - Analysis of GLIPR1L1 expression in transfected 293T cell line. Lane 1: MGC26856 transfected lysate (25.63 KDa). Lane 2: Non-transfected lysate.



Flow Cytometry: GLIPR1L1 Antibody [H00256710-B03P] - Analysis of negative control 293 cells (Black) and GLIPR1L1 expressing 293 cells (Green) using GLIPR1L1 purified mouse polyclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00256710-B03P

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97019-5mg	Mouse IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00256710-B03P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

