

Product Datasheet

BOD1P Antibody (1C9) - Azide and BSA Free H00284257-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00284257-M01

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00284257-M01



H00284257-M01

BOD1P Antibody (1C9) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1C9
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Novus Biologicals Mouse BOD1P Antibody (1C9) - Azide and BSA Free (H00284257-M01) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	284257
Gene Symbol	BOD1L2
Species	Human
Specificity/Sensitivity	This product is specific for Human FAM44C monoclonal antibody (M01), clone 1C9 [Gene ID: 284257].
Immunogen	FAM44C (AAH21740.1, 1 a.a. ~ 172 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MADGGGGGSGGAGPASTRASGGGGPINPASLPPGDPQLIAIVGQLKSRGLFD SFRRDCKADVDTKPAYQNLSQKADNFVSTHLDKQEWNPANDNQLHDGLRQ SVVQSGRSEAGVDRISSQVVDPKLNHIFRPQIEQIIHEFLVAQKEAAVPALPPEP EGQDPPAPSQDTS
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Sandwich ELISA
Recommended Dilutions	Western Blot, ELISA, Sandwich ELISA
Application Notes	Mouse monoclonal antibody raised against a full-length recombinant FAM44C.




Images

Western Blot: BOD1P Antibody (1C9) [H00284257-M01] - Analysis of FAM44C expression in transfected 293T cell line by FAM44C monoclonal antibody (M01), clone 1C9. Lane 1: FAM44C transfected lysate (Predicted MW: 18.1 KDa). Lane 2: Non-transfected lysate.



Sandwich ELISA: BOD1P Antibody (1C9) [H00284257-M01] - Detection limit for recombinant GST tagged FAM44C is 0.3 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00284257-M01

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00284257-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

