

Product Datasheet

WT1 Antibody (960522) - Azide and BSA Free MAB5729

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/MAB5729

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/MAB5729



MAB5729

WT1 Antibody (960522) - Azide and BSA Free

Product Information	
Unit Size	100 ug
Concentration	1.0 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	960522
Preservative	No Preservative
Isotype	IgG2b
Purity	Protein A or G purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Mouse WT1 Antibody (960522) - Azide and BSA Free (MAB5729) is a monoclonal antibody validated for use in Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	7490
Gene Symbol	WT1
Species	Human
Immunogen	E. coli-derived recombinant human WT1 Met127-Gly249 Accession # P19544
Notes	0

Product Application Details	
Applications	Flow Cytometry, CyTOF-ready, Immunocytochemistry
Recommended Dilutions	Flow Cytometry, Immunocytochemistry, CyTOF-ready
Application Notes	This antibody is Cytof ready.




Images

Immunocytochemistry: WT1 Antibody (960522) [MAB5729] - WT-1 was detected in immersion fixed K-562 cell line using MAB5729 at 8 ug/ml for 3 hours at room temperature. Cells were stained using NL007 (red) and counterstained with DAPI (blue). Specific staining was localized to cytoplasm & nuclei.



Flow Cytometry: WT1 Antibody (960522) [MAB5729] - An intracellular stain was performed on THP-1 cells with WT1 [960522] Antibody MAB5729AF647 (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 2.5 ug/mL for 30 minutes at room temperature. Both antibodies were conjugated to Alexa Fluor 647.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to MAB5729

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
H00007490-Q01-10ug	Recombinant Human WT1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/MAB5729

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

