

Product Datasheet

P2X7/P2RX7 Antibody (545205) [Unconjugated] - BSA Free MAB97901

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/MAB97901

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/MAB97901



MAB97901

P2X7/P2RX7 Antibody (545205) [Unconjugated] - BSA Free


Product Information	
Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	545205
Preservative	0.02% Sodium Azide
Isotype	IgG2a
Conjugate	Unconjugated
Purity	Protein A or G purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Mouse P2X7/P2RX7 Antibody (545205) [Unconjugated] - BSA Free (MAB97901) is a monoclonal antibody validated for use in IHC. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5027
Gene Symbol	P2RX7
Species	Human
Immunogen	E. coli-derived partial recombinant human P2X7/P2RX7 protein (amino acids 356-595)

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 5 ug/ml, Immunohistochemistry-Paraffin 5 ug/ml

Images

Immunohistochemistry-Paraffin: P2X7/P2RX7 Antibody (545205) [Unconjugated] [MAB97901] - P2X7/P2RX7 was detected in immersion fixed paraffin-embedded sections of human bladder using Mouse Anti-Human P2X7/P2RX7 Monoclonal Antibody (Catalog # MAB97901) at 5 ug/mL for 1 hour at room temperature followed by incubation with the Anti-Mouse IgG VisUCyte HRP Polymer Antibody (Catalog # VC001). Before incubation with the primary antibody, tissue was subjected to heat-induced epitope retrieval using Antigen Retrieval Reagent-Basic (Catalog # CTS013). Tissue was stained using DAB (brown) and counterstained with hematoxylin (blue). Specific staining was localized to cytoplasm and cell membrane of epithelial cells in convoluted tubules. View our protocol for IHC Staining with VisUCyte HRP Polymer Detection Reagents.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to MAB97901

MAB003	Mouse IgG2a Isotype Control (20102) [Unconjugated]
H00005027-P01-2ug	Recombinant Human P2X7/P2RX7 GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]
NBL1-14033	P2X7/P2RX7 Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/MAB97901

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

